



## Optical Alignment Deformation Spectroscopy

*Keith B. Neeves, David W.M. Marr, Kevin B. Roth, and Charles D. Eggleton*

**Summary:** A microfluidic device and method for high-throughput measurement of cell deformation


**Description:** Cell mechanical properties are a useful measure of phenotype that can be quantified by cell deformability. There is a lack of high-throughput methods to investigate the mechanical properties of large populations of individual cells. To address this need, researchers at Mines have developed optical alignment deformation spectroscopy (OADS), a technique where hydrodynamic interactions between individual cells are used to create deformation. In OADS, a linear optical trap is used to align two incoming cells in a microfluidic cross-flow geometry, allowing hydrodynamic forces to induce a collision between cells at the stagnation point. After the interaction, the cells leave the stagnation point and a new pair of cells enters the trap. A convenient model cell to characterize OADS is the human erythrocyte because of its well-known mechanical properties. Deformation data of erythrocytes is fit to a linear viscoelastic constitutive model (Voigt). The results show OADS has potential as an accurate, high-throughput, individual cell mechanical cytometer.

### Main Advantages of this Invention:

- High through-put assay measuring individual cells instead of population average
- Reagent and label free
- More accurate than current methods

### Potential Areas of Application:

- Disease screening
- High through-put blood testing method



Device schematic and view of cross-junction, with cells entering and exiting

**ID number:** 12012

**Publications:** K.B. Roth *et. al.*, *Lab Chip* 2013, 13, 1571. (Available upon request.)

**Intellectual Property Status:** US utility patent pending (application #13/770,875)

**Opportunity:** We are seeking an exclusive or non-exclusive licensee for the manufacturing, marketing, and sale of this technology.

### For more information contact:

William Vaughan, Director of Technology Transfer  
Colorado School of Mines, 1500 Illinois Street, Guggenheim Hall Suite 314, Golden, CO 80401  
Phone: 303-384-2555; e-mail: wvaughan@mines.edu